

Intrafix® SafeSet

Packaging information

1 Intrafix® SafeSet

2 Air Stop **3** Prime Stop **4** 15 µm **5** 20 ml **6** Latex-free **7** P

M.R.P. : ₹ 220.00 (Incl. of all taxes) 0617

8 YYYY-MM-DD 180 cm **13**

9 YYYY-MM-DD

10 LOT Batchno. **14** STERILE EO

11 REF 4063000 **15** DMV

12 40°C; 2 bar; 6,3 ml/1m GTIN: 04022495869517

16 **17** See instructions for use

18 Caution **19** For single use only **20** CE 0123 **21** Recycled **22** Do not use if package is damaged **23** PHT & DEHP-free **24** Non-pyrogenic **25** Conformity to Ukrainian standards **26** Conformity to Russian standards **27** Manufacture site

22 Do not use if package is damaged
23 PHT & DEHP-free
24 Non-pyrogenic
25 Conformity to Ukrainian standards
26 Conformity to Russian standards
27 Manufacture site

B. Braun Melsungen AG
34209 Melsungen, Germany
www.bbraun.com

B | BRAUN

1. Product name
2. AirStop
3. 20 drops per ml
4. Infusion filter 15 µm
5. PrimeStop
6. Latex-free
7. Pressure (no Gravity)
8. Expiration date
9. Date of manufacture
10. LOT/Batch number
11. Reference number
12. Maximum pressure
13. Tubing length
14. Sterilisation method
15. GTIN(GS1) code
16. Product picture
17. See instructions for use
18. Caution
19. For single use only
20. CE certification
21. Recycled
22. Do not use if package is damaged
23. PHT & DEHP-free
24. Non-pyrogenic
25. Conformity to Ukrainian standards
26. Conformity to Russian standards
27. Manufacture site